



Manufacturer	Alpaqua Engineering, LLC
Part #	A000020
Manufacturer's Website	http://www.alpaqua.com/Products/HeatingCooling/POGO.aspx
Description	Passive Heating/Cooling Block for 24 Tubes
Country of Origin	U.S.
Harmonized Code	8479909440 Parts for Industrial Robots NESOI
US Export Restrictions	No export license from the U.S. Department of Commerce is required, however, all restrictions governing export or re-export to denied countries, entities, and individuals apply. See http://export.gov for details.
Material	Aluminum, Acetal, Stainless Steel
Dimensions (LxWxH)	Product: 127.8 mm x 85.6 mm x 42.3 mm / 5.03" x 3.36" x 1.66" Shipping: 195 mm x 135 mm x 80 mm / 7.75" x 5.25" x 3.25"
Weight	Net Weight: 1000 g / 2.2 lb Shipping Weight: 1150 g / 2.5 lb
MSDS	Available upon request
RoHS/REACH	This product is compliant with RoHS/REACH regulations in effect on the date shown below.
Additional Comments	This product does not fall within the scope of products covered by EU directives governing the requirements for CE marking and therefore does not bear a CE mark. Alpaqua Products are intended to be used for research purposes only.

A000020 Product Passport.Docx - Last update: July 30, 2026