




Manufacturer	Alpaqua Engineering, LLC
Part #	A000280
Manufacturer's Website	https://www.alpaqua.com/product/24-position-sbs-tube-rack/
Description	24-Position Tube Rack (closed bottom)
Country of Origin	U.S.
Harmonized Code	9033.00.0000 Parts for Chemical Analysis Instruments
US Export Restrictions	No export license from the U.S. Department of Commerce is required, however, all restrictions governing export or re-export to denied countries, entities, and individuals apply. See http://export.gov for details.
Material	Acrylic, Aluminum, Stainless Steel
Dimensions (L x W x H)	Product: 127.8 mm x 85.6 mm x 67.25 mm / 5.03" x 3.36" x 2.65" Shipping: 195 mm x 135 mm x 105 mm / 7.75" x 5.25" x 4.13"
Weight	Net Weight: 235 g / 0.5 lb Shipping Weight: 360 g / 0.8 lb
MSDS	Available upon request
RoHS/REACH	This product is compliant with RoHS/REACH regulations in effect on the date shown below.
Additional Comments	For Research Use Only

A000280 Product Passport.Docx - Last update: July 30, 2026