



<b>Manufacturer</b>	Alpaqua Engineering, LLC
<b>Part #</b>	A000280
<b>Manufacturer's Website</b>	<a href="https://www.alpaqua.com/product/24-position-sbs-tube-rack/">https://www.alpaqua.com/product/24-position-sbs-tube-rack/</a>
<b>Description</b>	24-Position Tube Rack
<b>Country of Origin</b>	U.S.
<b>Harmonized Code</b>	9033.00.0000 Parts for Chemical Analysis Instruments
<b>US Export Restrictions</b>	No export license from the U.S. Department of Commerce is required, however, all restrictions governing export or re-export to denied countries, entities, and individuals apply. See <a href="http://export.gov">http://export.gov</a> for details.
<b>Material</b>	Acrylic, Aluminum, Stainless Steel
<b>Dimensions (W x D x H)</b>	6" x 4" x 4" 155 x 97 x 105 mm
<b>Weight</b>	0.52 lbs., 235 g
<b>MSDS</b>	Available upon request
<b>RoHS/REACH</b>	This product is compliant with RoHS/REACH regulations in effect on the date shown below.
<b>Additional Comments</b>	

A000280 Product Passport.Docx - Last update: January 9, 2023